

medin **SINDI**[®] nCPAP driver

Comprehensive respiratory support for pretermes and newborns from the first minute of life.



Mode nCPAP
with Medijet[®]



- Variable Flow
- Flow controlled up to 16 l
- Modified benveniste valve
- Pressure and O₂ measurement with alarm

Mode HHHFNC
Heated humidified
high flow nasal canula



- O₂ measurement with alarm
- Change to Medijet[®] by using same circuit

Giving future...

medin[®]
medical innovations

Circuits for Medijet® and Nasal Canulas

REF	Name	Inspiration/ expiration	Heated	Accessories	Single use/ reuse	Packing unit
5802	Patient circuit	Inspiration	yes	chamber included	Single use	15

Medijet®




REF	Name	Single use/ reuse
1000	Medijet® (20pieces/unit)	Single use
1000-10	Medijet® Starter Set (for 10 patients)	Single use
1000-20	Medijet® Starter Set (for 20 patients)	Single use

For detailed information concerning the Starter Sets please look at our product catalogue or contact our customer support.



Prongs

Size	x-small	small	medium	large	x-large	medium wide	large wide
REF	1200-01	1200-21	1200-02	1200-22	1200-03	1200-32	1200-33
Nostril Ø [mm]	3	3.5	4.1	4.75	5.5	3.7	5
Packing unit	10	10	10	10	10	10	10



Masks

Size	small	medium	large	x-large
REF	1200-04	1200-05	1200-06	1200-07
Mask width nose [mm]	18.5	21.0	23.0	26.0
Packing unit	10	10	10	5



Bonnets - Disposable

Size	XXS	XS	S	M	L	XL	XXL	XXXL
REF	1213-10	1214-10	1215-10	1216-10	1217-10	1218-10	1219-10	1220-10
Head circumference [cm]	17 - 19	19 - 21	21 - 23	23 - 25.5	25.5 - 28	28 - 30	30 - 33	33 - 36
Packing unit	10	10	10	10	10	10	10	10

Comfort Soft Plus



REF	QTY	Name	Size	Flow rate
1390	3/size	Set with three sizes (premature, neonatal, paediatric)	all three	depends on size

www.medingmbh.com



medin Medical Innovations GmbH
 Adam- Geisler- Str. 1
 82140 Olching, Germany
 Phone: +49 8142 448 460
 info@medingmbh.com



REV01

Subject to change without notice.